

Advantage® HD Polymerase Mix Protocol-at-a-Glance

(PT3921-2)

Advantage HD Polymerase (Cat. No. 639241) is a novel DNA polymerase that offers both high accuracy and high efficiency for PCR amplification of DNA. Its unmatched performance is achieved by its superior proofreading function, due to a robust 3'→5' exonuclease activity. Additionally, high specificity is achieved by inclusion of a Hot Start antibody that inhibits polymerase activity at ambient temperatures, preventing false initiation events during reaction assembly due to mispriming or primer digestion. When used with Clontech's optimized reaction buffer, Advantage HD Polymerase achieves the high fidelity, high sensitivity, and high specificity required for critical applications such as gene cloning and amplification of cDNAs for library construction.

The key features of this enzyme mix are high accuracy, with only 12 errors per 250,000 bp, rapid annealing, and improved specificity due to increased priming efficiency. It also exhibits higher amplification efficiency than *Taq* Polymerase and excellent performance even on GC-rich templates. This allows the use of a single PCR cycling protocol for products of varying sizes, including up to 8.5 kb on human genomic DNA, 10 kb on *E. coli* genomic DNA, and 22 kb on λ DNA.

The protocol for using Advantage HD Polymerase Mix is outlined below.

PCR Protocol

1. Thaw all reagents on ice.
2. Assemble the PCR reactions at room temperature, by adding each of the components listed in Table I to tubes or plates.

TABLE I. COMPONENTS OF PCR REACTIONS	
Reagents	Volume/Final Concentration per Reaction
Sterile deionized H ₂ O	up to 25 μ l
5X Advantage HD Buffer (Mg ²⁺) ^a	5 μ l
dNTP Mixture (2.5 mM each)	2 μ l
Primer 1	0.2–0.3 μ M
Primer 2	0.2–0.3 μ M
Advantage HD Polymerase (2.5 units/ μ l) ^b	0.25 μ l
Template	200 ng

^a Each vial contains 1 ml of 5X Advantage HD Buffer (5 mM MgCl₂).

^b The enzyme is supplied in a storage buffer containing 50 mM Tris-HCl (pH 8.2 at 4°C), 100 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 0.1% Tween®20, 0.1% Nonidet P-40, and 50% glycerol.

3. Mix well and spin down briefly to collect all the liquid at the bottom of the wells.



Clontech

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

(PR692030; published 20 September 2006)

4. Use one of the following PCR Protocols.

Note: Clontech recommends first trying the Three-Step PCR Method with Advantage HD Polymerase.

Three-Step PCR Method		Two-Step PCR Method	
30 cycles:		30 cycles:	
98°C	10 sec ^a	98°C	10 sec
55°C ^b	5 sec or 15 sec ^c	68°C	1 min/kb
72°C	1 min/kb		

^a Denaturation Conditions: Denature at 98°C for 5–10 sec. If a lower denaturation temperature of 94°C is used, increase time to 10–15 sec.

^b Annealing Temperature: Anneal at 55°C initially. Annealing temperature optimization may be required.

^c Annealing Time: Annealing time is dependent on primer length and on T_m values. For primers <25mers, anneal for 5 sec. For primers ≥25mers, calculate T_m values using the following formula, and choose the appropriate annealing time:

$$T_m \text{ value (°C)} = 2(\text{NA} + \text{NT}) + 4(\text{NG} + \text{NC}) - 5$$

When $T_m \geq 55^\circ\text{C}$, anneal for 5 sec.

When $T_m < 55^\circ\text{C}$, anneal for 15 sec.

Important Note: Advantage HD Polymerase has an extremely high priming efficiency. Annealing times longer than 5–15 sec may cause smearing of PCR products. Please try the Two-Step PCR Method when smeared DNA products are observed with the Three-Step PCR Method, or when using primers with T_m values ≥ 70°C.

PCR products obtained using Advantage HD Polymerase cannot be used directly for TA cloning. The termini of PCR products obtained with Advantage HD Polymerase are blunt-ended due to this enzyme's 3'-5' exonuclease activity. Advantage HD Polymerase products should be blunt-end cloned into appropriate vectors. When performing blunt-end cloning using a dephosphorylated vector, PCR products should also be phosphorylated, or PCR primers having phosphoric acid residues at their 5' termini must be used.

General Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

NOTICE TO PURCHASER: LIMITED LICENSE

Use of this product is covered

6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim (such as method claims in U.S. Patent Nos. 5,994,056 and 6,171,785) and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby conveyed by the purchase of this product expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

The Hot Start antibody is licensed under U.S. Patent No. 5,338,671.

Tween® is a trademark of ICI Americas, Inc.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc.

Clontech is a Takara Bio Company. ©2006

Produced by Takara Biotechnology (Dalian) Co., Ltd.