
RNA to cDNA EcoDry™ Premix (Double Primed)**Catalog No.**

639570

Amount

8 rxns

Lot Number

Specified on product label.

Description

RNA to cDNA EcoDry Premix (Double Primed) is a convenient dry master mix that allows efficient and accurate first-strand cDNA synthesis. Each well contains a mix of SMART™ MMLV Reverse Transcriptase, dNTPs, buffer, and both random hexamer and oligo(dT)₁₈ primers. Simply reconstitute the master mix by adding PCR-grade water along with your RNA. EcoDry Premixes are packaged as individual 8-well tube strips with optically clear sealing caps.

Package Contents

- 1 x 8-tube strip RNA to cDNA EcoDry Premix (Double Primed)
- 1 Optically Clear PCR Cap Strip (8 caps/strip)

Storage Conditions

- Always store any unused product at 20–22°C, either (1) back in the original foil pouch, re-sealed with desiccant; or (2) in a desiccator.

Shelf Life

- 1 year from date of receipt under proper storage conditions.

Shipping Conditions

- Room temperature

Product User Manuals

User manuals for Clontech® products are available for download at www.clontech.com/manuals.

The following user manuals apply to this product:

- RNA to cDNA EcoDry Premix Protocol-at-a-Glance (PT5153-2)

Quality Control Data**Raw Material Quality Control**

A sample aliquot of raw material SMART MMLV Reverse Transcriptase was tested to ensure optimal performance by first strand cDNA synthesis and dNTP incorporation assays. The purified enzyme was shown to be free of detectable levels of endonuclease, exonuclease, RNase and DNA contamination.

Functional Quality Control

For each lot, three tubes from a single strip were used to perform first-strand cDNA synthesis. In this assay, template RNA (a synthetic poly A⁺ RNA ladder ranging in size from 0.5 to 10 kb) was added to the RNA to cDNA EcoDry Premix and incubated for 1 hour at 42°C. The single-stranded cDNA products were separated on an alkaline-agarose gel, stained with a single-strand DNA binding dye, and visualized with UV light; the expected pattern was observed.

Clontech Laboratories, Inc.

A Takara Bio Company

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: tech@clontech.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (050113)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.543.6116

RNA to cDNA EcoDry™ Premix (Double Primed) (Sample-Not for Sale)

CATALOG NO.

639570

NOTICE TO PURCHASER:

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements listed below and described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

TRADEMARKS:

Clontech, the Clontech logo, EcoDry, and SMART are trademarks of Clontech Laboratories, Inc.

All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2013 Clontech Laboratories, Inc. This document has been reviewed and approved by the Clontech Quality Assurance Department.

Clontech Laboratories, Inc.

A Takara Bio Company

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: tech@clontech.com

5/1/2013

United States/Canada

Asia Pacific

Europe

Japan

800.662.2566

+1.650.919.7300

+33.(0)1.3904.6880

+81.(0)77.543.6116