


## pTRE3G-BI-mCherry

### Catalog No(s).

631333 (Not sold separately).

Sold as part of 631338, 631344, & 631341.



3' end of TRE3G-BI promoter

					NotI	NaeI	
		BglII					
	ApaI		BamHI		EagI		EcoRV
450	AAGTCGACAC	CGGGGCCAG	ATCTCCGCGG	GGATCCATCG	ATCCGCGGCC	GCCGGCGATA	TCT
	TTCAGCTGTG	GCCCGGGTC	TAGAGGCGCC	CCTAGGTAGC	TAGGCGCCGG	CGGCCGCTAT	AGA

### pTRE3G-BI-mCherry vector map and multiple cloning site.

#### Clontech Laboratories, Inc.

A Takara Bio Company

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: [tech@clontech.com](mailto:tech@clontech.com)

United States/Canada

800.662.2566

Last update: 040212

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.543.6116

